

Confused About Zirconia Choices? Try Our Monolithic Zirconias




Imagine



769 mPA

Lava Esthetic



800 mPA

ZR PLUS



1,000–1,100 mPA

ZR



1,100 mPA

ZR-V



1,100 mPA

Anterior Choices

- Lava Esthetic typically better for Vita shades A1, A2, B1 and B2
- Imagine typically better for Vita shades A3, A4, Cs and Ds
- Up to 3-unit bridges; not for Bruxers

All-Around Choice

- More translucent than ZR or ZR-V with added strength



Posterior Choices

- ZR has strength with added incisal enhancer
- ZR-V has solid color throughout


Zirconia Preparation Guidelines

Over

Anterior Crown





BUTT MARGIN




CHAMFER

Posterior Crown



BUTT MARGIN



CHAMFER

