



Locator Flow Chart: New Denture

*** Pick up in the morning counts as first day in lab; pick up in the afternoon counts as next day in lab*

First Appointment

- Remove healing abutment – measure depth of tissue, from the top of the prosthetic head (buccal, mesial, distal, lingual)
- Use this measurement to order abutments and impression copings
- Choose locator abutment that is the height of your greatest depth
- Replace healing abutments
- Take impression for custom tray (3 days in lab)

Second Appointment

- Remove healing abutment
- Place locator abutments using a locator tool
- Take impression, using locator impression copings
- Replace healing abutments, paying special attention to abutment location if abutments are different lengths
- Send impression and all retentive males to lab, requesting base plate and wax rims on Rx (2 days in lab)

Third Appointment

- Take bite registration and return for set-up (7 days in lab)

Fourth Appointment

- Try-in wax-up and verify bite
- Return for proses (4 days in lab)

Fifth Appointment

- Insert denture – torque locator abutments to 30Ncm

*** Denture will be returned with a 3.0 lb. pink retentive male unless specified. The clear 5.0 lb. and blue 1.5 lb. are also included.*

Please contact us with any questions at (800) 236-3859