

# Lucitone Denture

## PATIENT CARE & CLEANING INSTRUCTIONS



- Remove your dentures every night. Do not wear them to bed. This allows the gum tissue underneath to rest.
- Brush your denture daily with a soft bristle toothbrush using a mild soap like Ivory. Then rinse well and STORE DRY in a container overnight.
- You can also use popular denture cleansers like Polident and Efferdent to clean your denture, but DO NOT SOAK OVERNIGHT in the cleaning solution or in water.
- Do not use HOT OR BOILING WATER as this may warp the denture.

➤ **QUESTIONS?** CONTACT YOUR DENTAL OFFICE