



Over-Denture Bar Flow Chart

*** Pick up in the morning counts as first day in lab; pick up in the afternoon counts as next day in lab*

First Appointment

- Take alginate impression
- Submit Rx requesting custom tray (3 days in lab)

Second Appointment

- Remove healing abutments – tools can be rented if needed
- Place impression campings, ensuring that the screw finger is tightened appropriately
- Take open tray impression for master cast
- Replace healing abutments
- Submit Rx requesting base plate and wax rims (3 days in lab)

Third Appointment

- Remove healing abutments
- Take bite registration
- Put verification jig in place and lute together
- Replace healing abutments
- Submit Rx requesting bar and set-up (3 weeks in lab)

Fourth Appointment

- Try-in bar for passive fit
- Verify set-up
- Return for processing (4 days in lab)

Fifth Appointment

- Seat bar and torque to proper specifications – tools can be rented if needed
- Insert denture

Please contact us with any questions at (800) 236-3859