







Traditional vs Partial Denture Flow Chart

** Pick up in the morning counts as first day in lab; pick up in the afternoon counts as next day in lab

TRADITIONAL DENTURE

First Appointment

 Take preliminary impressions to fabricate custom trays (3 days in lab)

Second Appointment

 Take master impression to fabricate baseplate and wax rim (2 days in lab)

Third Appointment

 Try in baseplate and wax rim and take bite registration; return for set-up (7 days in lab)

Fourth Appointment

 Try-in/tooth try in and return for processing (4 days in lab)

Fifth Appointment

• Insert full denture

PRINT DENTURE

First Appointment

 Take preliminary impressions to fabricate baseplate and wax rim (2 days in lab)

Second Appointment

 Take bite registration using baseplate and wax rim and then take reline wash final impression inside the baseplate, using it as your custom tray; return for printed try in (3 days in lab)

Third Appointment

 Try-in and return with any adjustment notes for final printing (2 days in lab if no changes, 3 days if adjustments needed)

Fourth Appointment

· Insert full denture